

TECHNICAL DATA SHEET

PARAFFINS | SOLID

Pharmaceutical-Cosmetic | Technical Grades

Paraffins	Congealing Point	Oil Content	Needle Penetration	Kinematic Viscosity
	ISO 2207 at 15°C	ISO 2908 % w/w	ASTM D1321 Mm x 0.1 @25°C	ASTM D445 Mm ² /s @100°C

Fully Refined Paraffin Waxes

Wax 50/52 DAB	51	0.4	18	3.2
Wax 52/54 DAB	53	0.4	18	3.5
Wax 54/56 DAB	55	0.4	18	3.9
Wax 56/58	57	0.4	18	4.2
Wax 58/60 DAB	59	0.5	18	4.5
Wax 58/60 D DAB	58	0.4	14	4.2
Wax 60/62 STD DAB	61	0.5	18	4.9
Wax 60/62 321 DAB	61	0.4	15	4.7
Wax 62/64 P 476	63	0.7	20	5.9
Wax 64/66	66	0.5	19	6.2

Drawing Paraffines

Wax 54/56 P 472	55	0.7	22	4.2
Wax 56/58 P 527	57	0.5	18	4.0

Specialties

Wax 45/47	46	0.7	45	2.7
Wax 48/50	49	0.5	25	3.1
Wax 48/50 405 DAB	49	1.3	50	3.1
Wax 50/52 405 DAB	51	1.3	50	3.4
Wax 52/54 405 DAB	53	1.0	45	3.6
Wax 44/48 Z	46	2.0	50	3.0
Wax 48/50 Z	50	1.3	50	3.2
Wax 50/52 Z	51	0.3	40	3.2
Wax 52/54 Z	53	1.3	40	3.9
Wax P 144	42	1.0	60	2.45
Gepaf 6062	66	1.0	22	7.4
Ersolit 336 G	70	2.8	50	11.2
Wax P	72	2.5	50	17.0
Microfin 74	72	1.0	20	10.6

Data are subject to usual tolerance (mm²/s = cSt)